


PLAN D'ENSEMBLE



REP.	QTE.	NÁDE PIECE	DESCRIPTION
1	1	SF45042	INOX ELECTROPOLI
2	1	SF45043	INOX ELECTROPOLI
3	1	SF44675	SILICONE BLEU
4	1	SF44096-2	BLEU CLAIR
5	1	PA40005-5	SILICONE JAUNE
6	1	PA40005-3	SILICONE ROUGE
7	5	PA40006-12	BLEU FONCE
8	1	CS0111	VIS TCF DIN 84 M2X3 A2
9	2	CS0991	RIVET INOX 2,4x6
10	1	ML929	V.2
11	1	ML930	V.3

04	27/06/23			ENDO STAR 5 NICHROMINOX
03	27/09/19			
NÁ	Date	Visa	Ctrl	
 NICHROMINOX 18, rue des Frères Lumières F-69720 ST-BONNET-DE-MURE Tel : 04 78 74 04 15 Fax : 04 78 01 61 80				Dessine : DB
				PN182160
				Page : 1/1