




REP.	QTE.	NÂDE PIECE	DESCRIPTION
1	1	SF42715	INOX ELECTROPOLI
2	1	SF42716	INOX ELECTROPOLI
3	1	SF44368	BLEU FONCE
4	1	SF44096-2	BLEU CLAIR
5	1	PA40005-3	SILICONE ROUGE
6	1	PA40005-4	SILICONE VERT
7	1	PA40005-5	SILICONE JAUNE
8	12	PA40006-12	BLEU FONCE
9	1	CS0111	VIS TCF DIN 84 M2X3 A2
10	2	CS0991	RIVET INOX 2,4x6
11	1	ML571	V.2
12	1	ML572	V.4

05	29/06/23			ENDO STAR 12 NICHROMINOX	
04	27/09/19				
NÂ	Date	Visa	Ctrl	Dessine : RF	Page : 1/1
NICHROMINOX 18, rue des Frères Lumière F-69720 ST-BONNET-DE-MURE Tel : 04 78 74 04 15 Fax : 04 78 01 61 80				PN182180	