




REP.	QTE.	NÁDE PIECE	DESCRIPTION
1	1	SF42713	INOX ELECTROPOLI
2	1	SF42714	INOX ELECTROPOLI
3	1	SF44367-2	SILICONE BLEU
4	1	SF44096-2	BLEU CLAIR
5	1	PA40005-3	SILICONE ROUGE
6	1	PA40005-4	SILICONE VERT
7	1	PA40005-5	SILICONE JAUNE
8	8	PA40006-12	BLEU FONCE
9	1	CS0111	VIS TCF DIN 84 M2X3 A2
10	2	CS0991	RIVET INOX 2,4x6
11	1	ML569	V.2
12	1	ML570	V.4

04	20/06/23			ENDO STAR 8 NICHROMINOX		
03	27/09/19					
NÁ	Date	Visa	Ctrl			
 18, rue des Fr ^s Lum ⁱ - tel F-69720 ST-BONNET-DE-MURE Tel : 04 78 74 04 15 Fax : 04 78 01 61 80				Dessine : DB	PN182170	Page : 1/1