

ISOLATION GOWN

Medical Device, Class I

TECHNICAL DATA SHEET

English

FEATURES & BENEFITS

- Isolation gown intended to protect both health care professionals and patient from cross-contamination, and acts as a barrier against fluids and particles during routine medical or hygiene procedures
- Non-Sterile single-use
- Gown made of non-woven fabric with full back closure. Features neck ties, string belt and sleeves closed with an elastic
- Latex-free

PRODUCT DESCRIPTION

Product name: Isolation gown, green
Product type: Single-use
Material: Spunbonded Polypropylene
Size: 130 cm width x 140 cm length
Weight: 20 gsm
Packaging: Outer case of 10 polybags x 10 pcs. each
Origin: China



For professional use only.

INDUSTRIES

Medical / Dental

CERTIFICATION & STANDARD

Complies with the requirements of European Regulation 2017/745 relating to Medical Devices.

Manufactured under ISO 13485:2016 certified system.



ISOLATION GOWN

Medical Device, Class I

TECHNICAL DATA SHEET

English

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun.

Shelf life of the product is 5 years.

LOGISTICAL INFORMATION

Outer case specifications

REF	Size (mm)	EAN Code
977-0428	450 x 250 x 330	3701407400994

Inner case (bag) specifications

REF	Size (mm)	QTY/ Outer case	EAN Code
977-0428	420 x 250 x 32	10 polybags of 10 pcs.	3701407401007

